te@hnology for liviag

Provincial Respiratory Outreach Program • prop@technologyforliving.org • 1.866.326.1245 Technology for Independent Living Program • til@technologyforliving.org • 604.326.0175

Recommended Cleaning Schedule for Bilevel Users

Daily Weekly Monthly ☐ Empty water from the □ Clean humidifier pot, ☐ Clean or replace filter if humidifier pot every hose, mask and needed morning. headgear ☐ Check your supplies (Use warm soapy water. (filters, cushions for Hand wash. Rinse Well. ☐ Wipe mask with warm masks, nasal pillows, damp cloth every Do not use Antibacterial, etc.) add to your list for morning to remove oils, scented soaps) your next order. *It is etc. ☐ Inspect hoses & masks advisable to keep a for wear and replace if back-up mask and hose ☐ Refill humidifier pot

with distilled water on hand. necessary before bed. *With proper care, face masks, nasal masks and nasal pillows will maintain their performance and can last up to and Suction Equipment (if applicable to you) Daily Weekly Monthly □ Clean suction collection ☐ Inspect filter, replace as □ Empty suction bottle with hot soapy needed collection bottle and water, rinse well rinse ☐ Flush yankauer and tubing with water after each use to keep tubing patent ☐ Ensure your suction unit is plugged in and charging

beyond one year.

For more information: www.technologyforliving.org