## te¢hnol⇔gy for livig

Provincial Respiratory Outreach Program • prop@technologyforliving.org • 1.866.326.1245

Technology for Independent Living Program • til@technologyforliving.org • 604.326.0175

## **Recommended Cleaning Schedule for Ventilator Users**

Daily	Weekly		Monthly
Empty water from the humidifier pot every morning. Let humidifier pot	Clean ventilator circuit  Clean humidifier pot		Take stock of respiratory supplies, create a reorder list and maintain a three
Empty water trap in the	Clean suction collection bottle	_	month supply  Rotate bedside ventilator with back-up ventilator if
morning  Refill humidifier pot with	Clean suction catheters		applicable
Distilled water before bedtime	Check back-up ventilator for battery level, turn On to verify function		Ensure settings on all ventilator are up-to-date
Change inner cannula and clean dirty inner cannula	Check ventilator filters		Replace ventilator filters as needed
Empty suction machine collection bottle			
Change suction catheter			
Check ventilator alarms			
Check function and location of Ambu bag			
Plug in portable suction machine to recharge battery overnight			
Restock and review emergency/back-up supplies bag			

## For more information: www.technologyforliving.org